

22502

MedLife Exam Table

Features:

- *Seamless premium vinyl upholstery*
Fully welded, powder-coated, all-steel frame
- *Pneumatic backrest and sliding leg extension*
- *Pull-out stirrups that folds/slides away when not in use*
- *Non-Slip sliding footstep and adjustable feet for leveling*
- *Drawers are accessible from both sides of the table*
- *Table paper holder*
- *Weight capacity: 500 lbs*

Specifications:

Table:

- **Length: 58.5" x Width: 27.5" x Height: 34"**
- **Length w/ foot extended: 75"**

Step Stool:

- **Width: 18" x Length 16.5" x Height 9.5"**

Drawers:



- **Width: 16.8" x Length: 20" x Height: 4"**

Additional Measurements:

- **Head section: 31.5"**
- **Seat Section: 39.5"**



Colour Options:

- 22502: Black 
- 22503: Mocha 
- 22504: Blue 