

C3A/C6A

Video Colposcope




C3A/C6A

Video Colposcope



Patient management software



With its bright LED cold lighting system, high-resolution camera and electronic green filter, C3A/C6A Video Colposcope realizes full-screen display of precise cervical examination, offering crucial information for clinical decision making to the gynecologist.

> Features

- LED cold lighting system reproduces an original view of tissue
- Advanced and fast auto-focusing system ensures clear image
- Remote capture control and user-friendly keyboard layout for easy operation
- Electronic green filter enhances vascular image display without compromise of illumination
- Acetic acid and iodine reaction test timer to capture the right and precise image
- Real-time display of magnification
- Powerful image display and data management software

R-way®

The R-way® evaluation system reveals the variation of “color” on the cervix to interpret the performance of cervical diseases.

CIN2+

Saline application –
‘R for Red’
To observe the area with rich blood circulation, this reflects in red.

5% acetic acid-‘W for White’
To reveal increased nuclear activity, this reflects in white.

Green filter- ‘A for Abnormal’
To better observe abnormal vascular features, this defines as A.

Lugol’s iodine solution-
‘Y for Yellow’
To observe fields of no iodine uptake, this reflects in yellow.

* R-way® is not available for sale in the USA



Capture the right and precise image

Analyze colposcopic findings

Export report in PDF format





C3A/C6A
Video Colposcope

About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- *In-Vitro* Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



Solic Medical

Ph: 416-740-8885
info@solicmedical.com




EDAN

A world of potential