

# S573/EWS

## Phlebotomy Equipment

### Features:

#### S573-Phlebotomy Chair:

- *Flip arm positions in place when in use, then flips out of the way for easy exit*
- *Adjustable height and depth flip arm*
- *11/4" square, heavy duty, all-welded, tubular steel, frame*
- *Neutral Gray, powder-coated frame*
- *Black contoured padded seat, arched backrest, flip arm and armrests*
- *Adjustable foot levelers*
- *400 lbs. load capacity under normal use*
- *Easy-clean, durable vinyl*
- *Seat: Height: 21" x Width: 19" x Depth: 16.5"*

### Shipping Dimensions:

- **Height: 19" x Width: 19.5" x Length: 27.5"**
- **Weight: 46 lbs.**

#### Elbow Wedge:


- *Made of anti-microbial vinyl and durable high-density foam*
- *Easy-to-clean and waterproof*
- *Keeps elbow straight and elevated*

### Size Options:

- **EWS: 3.5"x 10"x 2"**
- **EWM: 4.25"x 10"x 2"**
- **EWL: 5.20"x 10"x 2"**



S573



EWS