

SonoTrax Vascular

Ultrasonic Pocket Doppler



Solic Medical
Passion to lead



EDAN

A world of potential

SonoTrax Vascular Ultrasonic Pocket Doppler

SonoTrax Vascular has been widely adopted as standard equipment for the assessment of peripheral vascular disease. The versatile, compact and portable SonoTrax Vascular allows fast and simple but effective assessment of both arterial and venous disease.

Features

- Easy-to-use
- Compact design
- Support 1.5 V AA battery and Rechargeable Ni-MH Battery
- Accurate blood flow detection with high-fidelity sound
- Support earphone output
- Support 4 MHz, 5 MHz, 8 MHz waterproof probes
- Interchangeable waterproof probe (IPX8)

Specifications

- Audio Output Power: 1 W
- Size: 34 (D) x 89 (W) x 141 (H) mm
- Weight: about 300 g (including the battery)



IPX8 Waterproof Standard



Hi-Fi Speaker



Adjustable Volume



Support Earphone Output



Solic Medical

Ph: 416-740-8885
info@solicmedical.com




EDAN

A world of potential